

ARIETTA 750 DeepInsight x



Value Your Core Ultrasound With "x"



IMAGE QUALITY

- DeepInsight Technology
- eFocusing PLUS
- Carving Imaging
- Detective Flow Imaging (DFI)

WORKFLOW

- iRecognize
- Auto Optimizer
- Protocol Assistant

TRANSDUCERS

- L52H Hockey Stick Transducer
- C421 Micro Convex Transducer
- L35 Linear Transducer

RADIOLOGY

- iATT
 - Steatosis Backscatter Index (SBSI)
- Shear Wave Measurement (SWM)
- Contrast Harmonic Imaging (CHI)

REAL-TIME VIRTUAL SONOGRAPHY

- Real-time Virtual Sonography (RVS)
 - DICOM Q/R Auto Search

- Active Registration
- Powered Anatomy
- 3D Sim-Navigator/E-field Simulator
- Body Motion Tracking
- Needle Tracking

C23RV/C23 Micro Convex Transducer

SURGERY

- L43LAP Laparoscopic Linear Transducer
- L51K/L43K Intraoperative transducer
- C42T Intraoperative transducer

WOMEN'S HEALTH

- Data Coordination with Mammography
- Real-time Tissue Elastography (RTE)
 - Auto Frame Selection (AFS)
 - Assist Strain Ratio (ASR)

CARDIOVASCULAR

- i2DTT
- iDGD
- iVascular

ERGONOMIC DESIGN

ARIETTA 750 DeepInsight x

Experience high-definition imaging powered by DeepInsight technology, along with advanced applications that elevate examination performance and streamline workflow. The ARIETTA 750 DeepInsight x integrates today's advanced features to support productivity in ultrasound examinations across highly specialized clinical fields.

IMAGE QUALITY

In addition to enhancing spatial and contrast resolution, the ARIETTA 750 DeepInsight x incorporates advanced imaging technologies to consistently deliver clear, high-quality images across the entire screen.



DeepInsight Technology / eFocusing PLUS / Carving Imaging
OFF



DeepInsight Technology / eFocusing PLUS / Carving Imaging
ON



DeepInsight Technology / eFocusing PLUS / Carving Imaging
OFF




DeepInsight Technology / eFocusing PLUS / Carving Imaging
ON


Detective Flow Imaging (DFI)

Balance the flow sensitivity and the suppression of body motion noise


Works with Wall Motion Reduction PLUS to selectively suppress body motion noise, enabling high-sensitivity and high-visibility imaging. Noise reduction is dynamically applied based on the motion level of each body part.



Wall Motion Reduction PLUS




Judge the scattering of signal intensity to reduce noise



Three Advanced Technologies for Superior Image Quality

1 DeepInsight Technology

Achieve high image visibility through rich gradation expression



2 eFocusing PLUS

Ensure uniform image extraction across the entire screen

3 Carving Imaging

Emphasize the tissue structure keeping the tissue expression and texture

* The technology was developed and designed using deep learning, one of AI technologies. The performance and accuracy of the system does not automatically change after activation.


WORKFLOW

iRecognize*1*2



Automatic Body Mark/Probe Mark placement


During abdominal examinations, the scanning plane is automatically recognized, and Body Mark/Probe Mark are set upon pressing the freeze button.



Auto Optimizer

Reduce the effort required for gain adjustment.

Auto Optimizer continuously adjusts gain in real time to maintain consistent image clarity—even when the scanning plane changes. This contributes to smoother, more efficient examinations.



Protocol Assistant

Streamline examinations by registering procedures and conditions in advance.

Button operations are significantly reduced, supporting a more efficient workflow. Additionally, the Guide View function allows reference images to be displayed, enhancing consistency and accuracy during scanning.




* 1 The technology was developed and designed using deep learning, one of AI technologies. The performance and accuracy of the system does not automatically change after activation.
* 2 The results of automated function need the confirmation of an examiner. The results can be modified manually as needed.

TRANSDUCERS

L52H High Frequency Hockey Stick Transducer

Flexible operation in any area

The versatile shape fits a variety of grips, and the high-frequency technology delivers high image quality in various areas. Supports DFI, enabling detailed blood flow display.



C421 Micro Convex Transducer

Reliable performance with an ergonomic, easy-to-scan design

DeepInsight technology and eFocusing PLUS deliver both exceptional image quality and intuitive operation. With Wide Scanning, the viewing angle is expanded without compromising frame rate performance.



L35 Linear Transducer

Enhancing Tissue Structure Visibility

The L35 linear transducer delivers high contrast resolution for clearer imaging of fine tissue structures. By utilizing single crystal technology and a matrix structure, it achieves improved penetration and enhanced tissue visibility.



iATT



Estimation of Hepatic Steatosis

The system displays ATT (Attenuation) index to estimate the degree of hepatic steatosis. To ensure measurement accuracy, multiple reflections from the body surface can be excluded, and ROI depth can be aligned with World Federation for Ultrasound in Medicine and Biology (WFUMB) guidelines, supporting consistent and reproducible results.

Steatosis Backscatter Index (SBSI)*

Diagnostic index that incorporates the backscatter coefficient and signal variation, in addition to ATT (Attenuation).

Adjustable ROI position



Steatosis Backscatter Index (SBSI)*

Shear Wave Measurement (SWM)

Tissue Stiffness Evaluation with SWM

SWM assesses tissue stiffness by calculating Vs, the shear wave propagation velocity. To ensure measurement reliability, it also provides a VsN indicator, supporting confident clinical decision-making.



Contrast Harmonic Imaging (CHI)

Support for Target Detection and Differential Diagnosis

Focus dependency is reduced by using eFocusing technology, therefore homogenous contrast image can be acquired from near to far field. Various contrast enhancement methods are supported, and they can be switched by single press of a button without interrupting the examination flow.



Low-MI



Pulse Inversion (PI)



Amplitude Modulation (AM)



CHI-eFLOW

* Index of liver steatosis is for reference purposes only and is not intended for diagnostic use.




REAL-TIME VIRTUAL SONOGRAPHY

Real-time Virtual Sonography (RVS)

RVS enables real-time fusion of ultrasound images with images from other modalities.

This allows direct comparison of lesions, leveraging the strengths of each imaging technique.



DICOM Q/R Auto Search Active Registraion*1*2



Automated workflows are expected to simplify image alignment

After RVS startup, volume data search and ultrasound-virtual image alignment are performed automatically, minimizing/reducing manual steps. A guidance display is also provided to support users throughout the process, ensuring a smooth and intuitive experience.

Active Regist. Capture Tips

1. Place probe on the body.
2. Select [Active Regist. Capture] button

3. Scan over 20 degrees in one direction.
4. Select [Active Regist. Capture] button or scan ends in 10 secs automatically.

RVSON

DICOM Q/R Auto Search
Data Import


Active Registration

Start

Powered Anatomy*1

Supporting Spatial Understanding


A region-classified reference image—generated through brightness analysis of CT data—is superimposed onto the virtual or ultrasound image.



3D Sim-Navigator/E-field Simulator

Estimate the puncture route based on target's shape

The spatial relationship between multiple electrodes (needles) and the target is visualized in 3D. Additionally, the electric field generated by radio waves—determined by electrode arrangement—can be superimposed onto CT or MRI images.



Body Motion Tracking

The omniTRAX Active Patient Tracker (manufactured by CIVCO) and the magnetic sensor keep the synchronization status of fused images even in case of misalignment by patient position change during CT/MR acquisition.

Needle Tracking

The VirtuTRAX Bracket (manufactured by CIVCO) and the magnetic sensor enable the tracking and display of the needle tip location during RFA. Additionally, the ARIETTA 750 DeepInsight x supports the function to correct needle tip displacement caused by needle bending and display it.

C23RV/C23 Micro Convex Transducer

Versatile support from routine to detailed examinations

The adoption of single crystal technology and an appropriate design --including a heat dissipation structure --is aimed to deliver the image quality comparable to our standard convex transducers and less stressful operation together.



* 1 The technology was developed and designed using deep learning, one of AI technologies. The performance and accuracy of the system does not automatically change after activation.
 * 2 The results of automated function need the confirmation of an examiner. The results can be modified manually as needed.

L43LAP Laparoscopic Linear Transducer

Flexible operability and support for puncture and ablation

Designed with a short-curvature mechanism inspired by endoscopic technology, the L43LAP is expected to offer superior operability for more efficient examinations. It is also equipped with two puncture notches to support precise puncture and ablation procedures.



L51K/L43K Intraoperative Transducer

Compact, User-Friendly Design for Surgical Support

Specifically designed for robotic-assisted surgeries in two different functional length. The flat contact surface delivers clear imaging - even on uneven organ surfaces. Also compatible with a wide range of laparoscopic forceps for stabler handling.



C42T Intraoperative Transducer

Combining High-Definition Imaging with Stable Scanning


Ergonomically designed to be comfortably held between the fingertips, the C42T enables steady and precise scanning during procedures. Equipped with advanced imaging technologies such as DeepInsight, it delivers clear, high-definition images across a wide range of clinical scenarios.



Real-time Tissue Elastography (RTE)


Real-time color display of tissue elasticity

In breast imaging, a series of operations - from selecting the appropriate measurement frame to initiating FLR (Fat Region Ratio) analysis - can be performed automatically, enabling easy and reproducible evaluations.



Data Coordination with Mammography

Mammography data can be viewed directly on the ultrasound system, reducing the need to move between devices during an examination. This integration helps maintain workflow continuity and is expected to streamline examinations.




CARDIOVASCULAR

i2DTT*1



Smooth workflow is expected to achieve faster measurement time

The operation flow is aligned with the examination process, enabling efficient and seamless procedures. The advanced i2DTT function allows simultaneous analysis of two-, three-, and four-chamber images, is expected to reduce overall measurement time.




iDGD*1*2




Support cardiovascular measurement with automation features and measurement

The system automatically measures E/e'-an index of left ventricular diastolic function-when using Dual Gate Doppler. Combined with the R-R Navigation function, it also automatically detects heartbeats suitable for measurement.



Measurement support functions*1


- Doppler Cursor Assist: Automatically detects the sample gate location to simplify Doppler cursor placement.
- Automatic Measurement: Performs EF (Ejection Fraction) measurement automatically-an essential index for evaluating cardiac function.
- Beat Mode: Automatically detects End-Diastole (ED) and End-Systole (ES) phases.
- R-R Navigation: Automatically identifies stable R-R intervals. Stability is defined by an interval ratio of approximately 1 between the preceding (RRp) and pre-preceding (RRpp) R-R intervals.



iVascular

Streamlined Doppler Measurement Support

During vascular examinations, settings such as the Color ROI position and cursor can be configured with a single click. This enables smooth transitions to subsequent Doppler measurements.



ERGONOMIC DESIGN

Appearance



10.4inch Touch Panel



23.8 inch LCD



Cable Assist Hook



Control Panel


* 1 The technology was developed and designed using deep learning, one of AI technologies. The performance and accuracy of the system does not automatically change after activation.
* 2 The results of automated function need the confirmation of an examiner. The results can be modified manually as needed.



ARIETTA 750 DeepInsight x

- External appearance and specifications are subject to change without notice.
- All products require the regulatory approval of the importing countries/regions.
- For details on their availability, contact our local representative.
- For proper use of the system, be sure to read the operating manual prior to placing it into service.
- FUJIFILM and FUJIFILM logo are registered trademarks or trademarks of FUJIFILM Corporation.
- The company names, product names and other names described in this brochure are trademarks or registered trademarks of FUJIFILM Corporation or its subsidiaries.
- ARIETTA 750 DeepInsight x is one of the ARIETTA 750 series.
- Fujifilm has been developing AI technologies that can be used for medical image diagnosis support, medical workflow support, and maintenance services for medical equipment, and is developing technologies that can be used in these areas under the brand name "REILI".
- The "REILI" marked features described in this brochure were developed and designed using machine learning or Deep Learning, which are parts of AI technology. The performance and accuracy of the system do not automatically change after implementation.

FUJIFILM

 **FUJIFILM Corporation**
26-30, Nishiazabu 2-chome, Minato-ku, Tokyo 106-8620, Japan
<https://global.fujifilm.com>

Ref. No. FXX-US-1030E (SK-24-10-X) ©2024 FUJIFILM Corporation